

聯合骨科器材股份有限公司

會議記錄單

頁次: 1 / 3

時 間：	2018 年 8 月 1 日 (星 期 三) 13 時 00 分 ~ 14 時 00 分				
地 點：	HC_12_2A、KS_3_4A、TP_B、V401				
主 席：	倪仲達		記 錄：倪仲達		
會議名稱：	Short neck type conformity Stem 專案啟動會議		權責單位：UR3		
出席人員：	何方元、簡惠蓮、何浩君、嶺榮娟、林成遠、鄭雅文、倪仲達、王坤祥、林清祺、詹世朗、陳怡親、羅正展、劉育良、陳文全、陳昕怡、呂御伸、余柏勳、梁志明、羅嘉慶				

項次	會議內容	負責人	預計日期	實際完成日	追蹤結果																																				
1	產品的需求為： Short neck type Conformity stem ✓ Short neck type (#1-#3) ✓ Collared type ✓ HA coated surface($155\mu\pm35\mu\text{m}$) ✓ Ti6Al4V material ✓ Neck length 28.5mm 	--																																							
2	相關的品號與尺寸分佈如下： <table><tr><th colspan="8">(30)</th></tr><tr><th>spec. #</th><th>品號</th><th>品名</th><th>stem length</th><th>offset</th><th>neck angle</th><th>neck length</th><th>distal width</th></tr><tr><td>1</td><td>PE 11101401</td><td>Conformity stem, short neck, collared, #1</td><td>115.0 mm</td><td>31.0mm</td><td rowspan="3">135°</td><td rowspan="3">28.5mm</td><td>7.5mm</td></tr><tr><td>2</td><td>PE 11101402</td><td>Conformity stem, short neck, collared, #2</td><td>119.5 mm</td><td>31.5mm</td><td>8.1 mm</td></tr><tr><td>3</td><td>PE 11101403</td><td>Conformity stem, short neck, collared, #3</td><td>124.0 mm</td><td>32.5mm</td><td>8.5 mm</td></tr></table>	(30)								spec. #	品號	品名	stem length	offset	neck angle	neck length	distal width	1	PE 11101401	Conformity stem, short neck, collared, #1	115.0 mm	31.0mm	135°	28.5mm	7.5mm	2	PE 11101402	Conformity stem, short neck, collared, #2	119.5 mm	31.5mm	8.1 mm	3	PE 11101403	Conformity stem, short neck, collared, #3	124.0 mm	32.5mm	8.5 mm	--			
(30)																																									
spec. #	品號	品名	stem length	offset	neck angle	neck length	distal width																																		
1	PE 11101401	Conformity stem, short neck, collared, #1	115.0 mm	31.0mm	135°	28.5mm	7.5mm																																		
2	PE 11101402	Conformity stem, short neck, collared, #2	119.5 mm	31.5mm			8.1 mm																																		
3	PE 11101403	Conformity stem, short neck, collared, #3	124.0 mm	32.5mm			8.5 mm																																		

聯合骨科器材股份有限公司

會議記錄單

頁次: 2 / 3

3	實質對等產品為 Corail stem(Depuy) 與MetaFix(Corin)。 實質對等產品 <table><tr><th>Product (manufacturer)</th><th colspan="2">Corail (DePuy)</th><th colspan="2">MetaFix (Corin)</th><th>Conformity stem (UOC)</th></tr><tr><td>Intended use</td><td colspan="2">Hip replacement</td><td colspan="2">Hip replacement</td><td>Hip replacement</td></tr><tr><td>Materials</td><td colspan="2">Ti-6Al-4V</td><td colspan="2">Ti-6Al-4V</td><td>Ti-6Al-4V</td></tr><tr><td>Surface finishing</td><td colspan="2">HA coating</td><td colspan="2">HA coating</td><td>HA coating</td></tr><tr><td>Collared</td><td colspan="2">V</td><td colspan="2">V</td><td>V</td></tr><tr><td>Collarless</td><td colspan="2">V</td><td colspan="2">X</td><td>X</td></tr><tr><td>Neck angle (°)</td><td>125</td><td>135</td><td>125</td><td>135</td><td>135</td></tr><tr><td>Stem length(mm)</td><td colspan="2">115-140</td><td>130-140</td><td>115-140</td><td>115-124</td></tr><tr><td>Neck length(mm)</td><td>35</td><td>32</td><td>30.5</td><td>31</td><td>28.5</td></tr><tr><td>Offset(mm)</td><td>34.8-39.6</td><td>33-34.2</td><td>38</td><td>33-33.5</td><td>31-32.5</td></tr><tr><td>FDA regulation compliance</td><td colspan="2">N/A</td><td colspan="2">K162942 K153381</td><td>---</td></tr><tr><td>CE regulation compliance</td><td colspan="2">CE 0459 CE 0086</td><td colspan="2"></td><td>---</td></tr></table>	Product (manufacturer)	Corail (DePuy)		MetaFix (Corin)		Conformity stem (UOC)	Intended use	Hip replacement		Hip replacement		Hip replacement	Materials	Ti-6Al-4V		Ti-6Al-4V		Ti-6Al-4V	Surface finishing	HA coating		HA coating		HA coating	Collared	V		V		V	Collarless	V		X		X	Neck angle (°)	125	135	125	135	135	Stem length(mm)	115-140		130-140	115-140	115-124	Neck length(mm)	35	32	30.5	31	28.5	Offset(mm)	34.8-39.6	33-34.2	38	33-33.5	31-32.5	FDA regulation compliance	N/A		K162942 K153381		---	CE regulation compliance	CE 0459 CE 0086				---	--			
Product (manufacturer)	Corail (DePuy)		MetaFix (Corin)		Conformity stem (UOC)																																																																								
Intended use	Hip replacement		Hip replacement		Hip replacement																																																																								
Materials	Ti-6Al-4V		Ti-6Al-4V		Ti-6Al-4V																																																																								
Surface finishing	HA coating		HA coating		HA coating																																																																								
Collared	V		V		V																																																																								
Collarless	V		X		X																																																																								
Neck angle (°)	125	135	125	135	135																																																																								
Stem length(mm)	115-140		130-140	115-140	115-124																																																																								
Neck length(mm)	35	32	30.5	31	28.5																																																																								
Offset(mm)	34.8-39.6	33-34.2	38	33-33.5	31-32.5																																																																								
FDA regulation compliance	N/A		K162942 K153381		---																																																																								
CE regulation compliance	CE 0459 CE 0086				---																																																																								
4	系統名:United Hip System 聯合髖關節系統 產品名: Conformity stem, short neck, collared, #1 (優康短頸有領股骨柄, #1)	--																																																																											
5	註冊預計申請地區: CE, FDA, TFDA, CFDA, PMDA ，CE認證與cementless type 一起送。	--																																																																											
6	建議應把conformity stem 系列的中英文品名統一，以英文表示方式為主，例如: Conformity stem, short neck, collared, #1 (優康短頸有領股骨柄, #1)	--																																																																											
7	標籤範例以下圖表示方式。 	--																																																																											
8	法規開時間與導量產時間需於後續確認是否需進行impingement 測試後提供	倪仲達	2018/8/24																																																																										
9	Stem 預計可以與Metal head以及ceramic head搭配	--																																																																											

聯合骨科器材股份有限公司

會議記錄單

頁次: 3 / 3

10	因與cementless type 一同評估並非worst case 目前暫無購買競品的需求。	--				
----	---	----	--	--	--	--